

Justification:	Function:	Electronic signature justification:	Signed: Date (yyyy-mm-dd) - Time (hh:mm):
Issued:	QA	Carolina Johansson - SEHRBJNC	2023-03-29 - 09:50
Reviewed:	QA	Abdallah Almashharawi - ABDALM	2023-03-29 - 09:59
Approved:	QA	Elin Andersson - ELIAND	2023-04-04 - 08:30
Released:	QA	Carolina Johansson - SEHRBJNC	2023-04-17 - 12:56

This document has been electronically signed by the persons above.

Released

Document No: 10000037228 Edition: 09 Release date: 2023-04-17

Product Information

Provox® BasePlate Adaptor


Product description:

The Provox BasePlate Adaptor ("adaptor") is an accessory product for rehabilitation after total laryngectomy. It allows attaching medical devices, (HME), with ISO 15mm standard connector to a tracheostoma by fitting it into a Provox Adhesive base plate, Provox LaryButton or Provox LaryTube. A typical example would be to attach an HME with built-in oxygen adapter (TrachPhone).

Provox BasePlate Adaptor facilitates the use of TrachPhone together with Provox Adhesives or Provox LaryTube.

Product Information

Document ID:	PF018-01-TechInfo	Edition:	09
Manufacturer:	Atos Medical AB Kraftgatan 8 SE-242 35 Hörby, Sweden		
Classification: (EU) 2017/745	Class I (Rule 1)		
Intended Use:	Provox BasePlate Adaptor is an accessory that allows attaching medical device, e.g. an HME, with an ISO 15 mm standard connector to a Provox attachment		
Use specifications:	Intended medical indication: Accessory product for patients after total laryngectomy Intended patient population: Male and female Typical average age: N/A. Cognitive ability, by a clinician judged as sufficient Manual dexterity: Unconscious patients must be constantly monitored. Not intended for patients with mechanical ventilation. Intended usage: Single patient use. Intended part of the body/type of tissue applied to or interacted with: Neck, (tracheostoma). Intended user profile: Patient, clinician, trained nurse. Intended conditions of use: Home use (normal daily environment without any hygienic or environmental restrictions regarding temperature, moisture etc.). Hospital use. Frequency of use: Continuous use. Replacement rate: Shall be changed after used for a maximum of 3 months.		
Contraindications:	Shall not be used for mechanical ventilation		
CE Mark:	Yes. Device is CE-marked		
GMDN code:	58705 (Tracheostoma protective filter)		
Sterilization:	Non-sterile		
Raw material:	Polyether ether ketone (PEEK)		
Latex information:	Not manufactured with natural rubber latex.		
Biological origin:	The device is not manufactured with materials derived from human or animal source.		
Handling and storage:	Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.		

Product Information

Waste handling and disposal: Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.

Hazardous components: None

Expiration date: 5 years after manufacturing.

Packaging: Provox BasePlate Adaptor is separately packed in a plastic bag of Low Density Polypropylene.
The products and instructions for use are packed in a cardboard box.

Devices under Basic UDI-DI: 7331791-HME-A-000-0003-F5

REF	Name	UDI-DI
7263	Provox BasePlate Adaptor	7331791001697

Atos Medical AB compatible products:

Range	BASIC UDI-DI
Provox LaryTube	7331791-LTU-0-000-0002-3E
Provox LaryButton	7331791-LTU-0-000-0000-38
Provox StabiliBase Provox XtraBase Provox StabiliBase OptiDerm Provox Flexiderm Provox Optiderm	7331791-ADH-0-000-0000-CQ
TrachPhone	7331791-HME-0-000-0006-XT
Freevent DualCare	7331791-HME-0-000-0005-XQ
Freevent XtraCare	7331791-HME-0-000-0004-XM